

BRAILLEX®

Live-Models

Fast work with notebook and smartphone or tablet



Presently, our daily life is defined by continuous availability.

Thanks to our three BRAILLEX Live-models, we have developed a line which meets the technical requirements of a blind person's everyday life. It is the perfect assistive device for every one in connection with a smartphone and/or a laptop. Discover our Live-family more closely.



FUNCTIONALITIES

Reading height/writing height
 Number of Braille cells
 Surface of cell caps
 Adjustable Braille dot firmness
 Routing keys per cell
 Type of keyboard
 4-step Easy Access Bar
 Modular operating modes of the Easy Access Bar
 Ports
 Operation mode
 optimized field of application

| | Live | Live+ | Live20 |
|--|-------------------------------|-------------------------------|---------------------|
| Reading height/writing height | 18 mm/23 mm | 18 mm/23 mm | 18 mm/23 mm |
| Number of Braille cells | 40 | 40 | 20 |
| Surface of cell caps | flat or concave | flat or concave | flat or concave |
| Adjustable Braille dot firmness | yes | yes | yes |
| Routing keys per cell | 1 | 1 | 1 |
| Type of keyboard | Braille | Qwerty | Braille |
| 4-step Easy Access Bar | yes | yes | yes |
| Modular operating modes of the Easy Access Bar | 3 | 3 | 3 |
| Ports | USB & Bluetooth | USB & Bluetooth | USB & Bluetooth |
| Operation mode | battery & cable | battery & cable | battery & cable |
| optimized field of application | smartphone, tablet & notebook | smartphone, tablet & notebook | smartphone & tablet |

Highlights:



- Simultaneous Operation of two devices by means of both ports
- Compatible with JAWS, SuperNova, NVDA, Apple VoiceOver and Google Talkback
- One-hand-operation thanks to the patented Easy Access Bar
- Pleasant reading experience and relaxed reading position

Technical data:



- (Live, Live+, Live20)
- Dimensions (w x l x h): 298 x 115 x 18/23 mm; 298 x 186 x 18/23 mm; 170 x 111 x 18/23 mm
- Weight: (incl. battery) 900 g; 1200 g; 547 g
- Number of Braille cells: 40; 40; 20
- Battery capacity: > 20 hours
- Compatibility: All common operating systems

Scope of delivery:

User's manual, case, power supply unit, USB connection cable

GO AND FIND OUT FOR YOURSELF!

Please feel free to contact:

Phone +49 2304 205 0

www.papenmeier-rehatechnik.de



BRAILLEX® Live, Live+, Live20

Photo Credits: @brechostudio - stock.adobe.com / @Jurgen Fälschle - stock.adobe.com / @F.H. Papenmeier GmbH & Co. KG R3755.617 08.2023
Changes and/or improvements of this product are subject to change without prior notice.