

# BRAILLEX<sup>®</sup>

## EL 40c

Mobile work place Equipment having  
a class of its own



**20 YEARS  
OF EXPERTISE**

**PAPENMEIER**



technologically leading - humanly close  
**PAPENMEIER**

RehaTechnik 

## BRAILLEX® EL 40c

Every Work place for the blind has different requirements. For flexible and mobile use in your daily job, we have developed our Braille display EL 40c. Connected to a notebook, it is the perfect assistive device for all blind users.



### FUNCTIONS

Reading height
Number of Braille cells
Surface of cell caps
Adjustable Braille dot firmness
Routing keys per cell
Freely selectable direct access keys
4-step Easy Access Bar
Modular operating modes of the Easy Access Bar
Port

### BRAILLEX® EL 40c

18 mm
40
flat or concave
yes
2
yes
yes
3
USB



#### Highlights:

- Accurate work: 40 high-quality Braille cells with integrated double routing keys to control cursor and mouse
- Increasing efficiency by means of freely selectable direct access keys
- One-hand-operation thanks to the patented Easy Access Bar
- Short familiarisation due to modular operation concept
- Pleasant reading experience and relaxed reading position
- Complete Braille display incl. Braille cells from German production



#### Technical data:

- Dimensions: (W x L x H) 298 x 92.5 x 18 mm
- Weight: 640 g
- Number of Braille cells: 40
- Compatibility: All common operating systems



#### Ships with the following components:


- BRAILLEX EL 40c
- Drivers for JAWS 32/64 bit, SuperNova, Apple VoiceOver, Linux
- User's manual
- Case
- USB connection cable

**GO AND FIND OUT FOR YOURSELF!**

**We would be glad to advise you at first-hand:**

Phone +49 2304 946 0

Internet [www.papenmeier-rehatechnik.de](http://www.papenmeier-rehatechnik.de)



BRAILLEX® EL 40c